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1. Laboratory biosafety guidance related to coronavirus disease (COVID-19) Interim guidance 13 May 2020 OMS
2. US Food and Drug Administration (FDA). Emergency Use Authorization (EUA) information, and list of all current EUAs. (<https://www.fda.gov/emergencypreparedness-and-response/mcm-legal-regulatory-andpolicy-framework/emergency-use-authorization>, accessed 3 May 2020).
3. Coronavirus disease (COVID-19) technical guidance: surveillance and case definitions. Geneva: World Health Organization; 2020 (<https://www.who.int/emergencies/diseases/novelcoronavirus-2019/technical-guidance/surveillance-andcase-definitions>, accessed 6 April 2020).
4. Laboratory biosafety manual: third edition. Geneva: World Health Organization; 2004 (<https://www.who.int/csr/resources/publications/biosafety/Biosafety7.pdf?ua=1>, accessed 6 April 2020).
5. CDC 2019-Novel Coronavirus (2019-nCoV) RealTime RT-PCR Diagnostic Panel. Atlanta: Centers for Disease Control and Prevention; 2020 (<https://www.fda.gov/media/134922/download>, accessed 22 April 2020)
6. Van Doremalen N, Bushmaker T, Morris DH, Holbrook MG, Gamble A, Williamson BN, et al. Aerosol and surface stability of SARS-CoV-2 as compared with SARS-CoV-1. *N Engl J Med*. 2020
7. Guidance on regulations for the transport of infectious substances 2019–2020. Geneva: World Health Organization; 2019 (WHO/WHE/CPI/2019.20; <https://www.who.int/ihr/publications/WHO-WHECPI-2019.20/en/>, accessed 6 April 2020).
- 8- Strengthening health security by implementing the International Health Regulations (2005). Biosafety video series. Geneva: World Health Organization. (<https://www.who.int/ihr/publications/biosafetyvideo-series/en/>, accessed 6 April 2020).
- 9- Banada PP, Sivasubramani SK, Blakemore R, Boehme C, Perkins MD, Fennelly K, Alland D. Containment of bioaerosol infection risk by the Xpert MTB/RIF assay and its applicability to point-of-care settings. *J Clin Microbiol*. 2010; 48(10):3551-7.
- 10- Tuberculosis laboratory biosafety manual. Geneva: World Health Organization; 2012 ([https://www.who.int/tb/publications/2012/tb\\_biosafety/en/](https://www.who.int/tb/publications/2012/tb_biosafety/en/), accessed 20 April 2020).
- 11- Chin A, Chu J, Perera M, Hui K, Yen HL, Chan M, et al. Stability of SARS-CoV-2 in different environmental conditions. *The Lancet Microbe* 2020. (<https://www.medrxiv.org/node/74521.externallinks.html>, accessed 7 April 2020).
- 12-. Kampf G, Todt D, Pfaender S, Steinmann E. Persistence of coronaviruses on inanimate surfaces and its inactivation with biocidal agents. *J Hosp Infect*. 2020; Feb 6. (<https://www.journalofhospitalinfection.com/action/showPdf?pii=S0195-6701%2820%2930046-3>, accessed 6 April 2020).

13. International Civil Aviation Organization (ICAO). Safety. Technical instructions for the safe transport of dangerous goods by air (Doc 9284) (<https://www.icao.int/safety/DangerousGoods/Pages/technical-instructions.aspx>, accessed 6 April 2020).
14. Managing epidemics: key facts about deadly diseases. Geneva: World Health Organization; 2018 (<https://apps.who.int/iris/handle/10665/272442>, accessed 6 April 2020).
15. How to handrub? With alcohol-based formulation. How to handwash? With soap and water. Geneva: World Health Organization; 2006 (<https://www.who.int/gpsc/tools/GPSC-HandRubWash.pdf>, accessed 6 April 2020).
- 16- Pan F, Xiao X, Guo J, et al. No evidence of severe acute respiratory syndrome-coronavirus 2 in semen of males recovering from coronavirus disease 2019. *Fertil Steril.* 2020;113(6):1135- 1139.
17. Diangeng Li, PhD; Meiling Jin, MD; Pengtao Bao, PhD; Weiguo Zhao, MD; Shixi Zhang, MD. Clinical Characteristics and Results of Semen Tests Among Men with Coronavirus Disease 2019, *JAMA Network Open.* 2020;3(5): e208292.
18. George Anifandis\*, Christina I. Messini, Alexandros Daponte, Ioannis E. Messinis COVID-19 and fertility: a virtual reality. *RBMO volume 00 issue 0 2020*
- 19.. Diangeng Li, PhD; Meiling Jin, MD; Pengtao Bao, PhD; Weiguo Zhao, MD; Shixi Zhang, MD. Weiguo Zhao, MD; Shixi Zhang, MD Clinical Characteristics and Results of Semen Tests Among Men With Coronavirus Disease 2019 *Research Letter | Infectious Diseases Jama Network Open* Mayo 7 (2020) Pag 1 a 3

**Escrito por: Julia Ariagno, Patricia Chenlo, Norma Pugliese, Melba Sardi, Herberto Repetto, Susana Curi y Gabriela Mendeluk.**

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